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We have been dealing
with air technology for more
than 30 years: specialization
has given good results.

6/A Via Natta. 43122 Parma. Italy
10°21'48" EST - 44°50'46" NORD

tel +39.0521.607604

fax +39.0521.607628 (Uff. vendite/ Sales Dept.)

fax +39.0521.607855 (Uff. acquisti/ Purchasing Dept.)

fax +39.0521.399966 (Uff. amministrativo/ Accounting Dept.)

www.cattani.it - e-mail: info@cattani.it

azienda con Sistema Qualità Certificato
UNI EN ISO 9001/2000 - UNI CEI EN ISO 13485

Company with Quality System certified according to
UNI EN ISO 9001/2000 - UNI CEI EN ISO 13485

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Respect for our environment
through efficient use of energy.
Sustainable progress.

Your suction system can also play a part in creating a better world.

For some time researchers have warned us that our present consumption of materials and energy is causing an excessive exploitation of our natural resources.

The less of these precious resources we are able to use, the greater is the chance that our ecosystem can regain its original capacity to provide them on a long-term and sustainable basis. On the other hand, the more of these resources we consume, the more our ecosystem's ability to supply them is diminished. Along with many of our choices as a company, our research and development has been guided by the belief that through consistent effort by us all our resources can be conserved for the benefit of future generations.

Micro-Smart is a suction system designed around the principles of eco-sustainability :

- reduction in materials used to produce the unit;
- the use of VSD (Variable Speed Drive) technology results in longer motor life, reduction in electrical consumption and operating costs;
- unique flexibility and new performance options to the clinic that were previously unavailable.

Micro-Smart's suction motor weighs only 7 kg, however thanks to VSD technology it outperforms suction motors of double its weight and dimensions. The VSD facilitates optimal utilisation of the motor and reduces operating costs.

The energy saving is most evident when the suction is underutilised; when the demand on the suction decreases the VSD reduces the power frequency to reduce speed. A reduction in motor speed results in a reduced energy absorption by the motor, providing a saving in electrical consumption.

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Micro-Smart is able to handle unexpected operational situations, without even stopping. Consistent suction is granted.

Thanks to the inverter control and computerised program Micro-Smart is able to handle a range of unexpected situations.

- 1** In difficult operational conditions, for example high temperatures, Micro-Smart continues to operate without any damage to the unit. It does this through automatically decreasing the vacuum level until the temperature returns to normal.
- 2** In the case of excessive fluids from the surgeries entering the Micro-Smart, an auto-protection system operates, momentarily opening a re-circulation valve, thus stopping the unit from being flooded. This also causes a momentary decrease in the speed of the suction motor (the system detects an increase in amperage to the centrifugal separator and reduces the speed of the suction motor), with the surgeries only noticing a minimal fluctuation if any at all, as the system goes back to its previous speed once the fluids are processed.
- 3** In cases of sudden voltage fluctuations and/or current overload within the system's parameters, the system continues to operate, thanks to Micro-Smart's computerised control.
- 4** Possible operating anomalies are shown on the display during operation. This 'self-diagnosis' greatly helps clinic staff and service technicians to identify and rectify an issue.

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70 mbar is recommended for the use of the saliva ejector to maintain a low level of fluids at the base of the mouth, without damaging delicate tissue

100 mbar is recommended for the use of the HVE spray interceptor which aspirates the stream of water and air from the syringe and rotating instruments, such as the sprays from scalers. Low pressure reduces discomfort and noise caused by the aspirated air to a minimum.

210 mbar recommended for the use of the surgical tip, which keeps the operating area clear and provides good visibility.

Adaptable performance.

Micro-Smart is designed to have three different preset pressure settings which have been programmed according to the feedback from our customers. Different settings however can be programmed by the operator, or the technician, in order to achieve the exact desired suction level, within 0 and 210 mbar.

By pressing the desired preset button on the touchpad control, the suction can be easily set according to the type of procedure to be performed.



The suction system can be installed outside of the surgery, whereas the touchpad control is mounted in the surgery (a telephone cable should be run between the Micro-Smart and the touchpad control).

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A display that communicates.

Through simple and easy operations on the keypad you can monitor the operating status of the suction system. The unit's software monitors a range of technical parameters, allowing precise analysis of any possible anomalies. The unit also has a memory function, which records a log of any errors.

Misure Kg = 23
 L = 380 mm
 P = 350 mm
 H = 600 mm
 Alimentazione: 230 Volt

Measurements Kg = 23
 W = 380 mm
 D = 350 mm
 H = 600 mm
 Supply: 230 Volt