




# SL Series Binocular Loupes



**Zumax Medical**  
*Beyond Vision*



## SLE Binocular Loupes



SLE Binocular Loupes



HL8200 LED HeadLight &  
SLE Binocular Loupes



Wooden Case

SLE series Binocular Loupes is designed for all possible medical needs.

- ◇ Multi-coated, high grade glass lenses produce crisp, clear image.
- ◇ The Loupes which employ eyepiece of very large view field diameter. The depth of focus is deep and object can be seen stereoscopically with actual distance feeling.
- ◇ The angle of the loupes can be adjusted between 0° ~40° with ease and fatigue is minimized even after prolonged use.
- ◇ 6mm up and down displacement. Pupil distance can be adjusted individually which could be suited for all kind of faces.
- ◇ Adapts to any Headband and Headlight of Zumax.
- ◇ Smart design, convenient for mobile diagnosis.
- ◇ Vertical distance: 6mm
- ◇ Adjustable viewing angle: 0° ~40°
- ◇ Adjustable pupil distance: 45mm/1.6"~80mm/3.1"


### Technical Specification SLE Binocular Loupes:

Magnification	Working distance (mm)	Field of View dia. (mm)
2.5X	340/420/460/500	80/100/110/120
3.0X	340/420/460/500	70/84/92/100
3.5X	340/420/460/500	60/72/78/84



### Zumax Medical Co.,Ltd.


36 Zhiying Street, Suzhou New District, Suzhou, Jiangsu 215129, China  
Tel: +86 512 666 50100 ext. 502/504/510 Fax: +86 512 6690 9655  
E-mail: sales@zumaxmedical.cn <http://www.zumaxmedical.com>



## SLH Binocular Loupes



SLH Binocular Loupes



SLH Binocular Loupes &  
HL8200 Headlight

SLH series Binocular Loupes is designed for all possible medical needs.

- ◇ Multi-coated, high grade glass lenses produce crisp, clear image.
- ◇ The Loupes which employ eyepiece of very large view field diameter. The depth of focus is deep and object can be seen stereoscopically with actual distance feeling.
- ◇ The angle of the loupes can be adjusted between 0° ~40° with ease and fatigue is minimized even after prolonged use.
- ◇ 6mm up and down displacement. Pupil distance can be adjusted individually which could be suited for all kind of faces.
- ◇ Adapts to any Headband and Headlight of Zumax.
- ◇ Smart design, convenient for mobile diagnosis.
- ◇ Vertical distance: 6mm
- ◇ Adjustable viewing angle: 0° ~40°
- ◇ Adjustable pupil distance: 45mm/1.6"~80mm/3.1"

### Technical Specification SLE Binocular Loupes:

Magnification	Working distance (mm)	Field of View dia. (mm)
2.5X	340/420/460/500	80/100/110/120
3.0X	340/420/460/500	70/84/92/100
3.5X	340/420/460/500	60/72/78/84



### Zumax Medical Co.,Ltd.

36 Zhiying Street, Suzhou New District, Suzhou, Jiangsu 215129, China  
 Tel: +86 512 666 50100 ext. 502/504/510 Fax: +86 512 6690 9655  
 E-mail: sales@zumaxmedical.cn <http://www.zumaxmedical.com>