



LOW SPEED HANDPIECES
ISO E-TYPE EXTERNAL WATER SPRAY

EX Series

versatile and affordable



NAKANISHI INC.
www.nsk-inc.com
700 Shimohinata Kanuma-shi
Tochigi 322-8666,
Japan
TEL : +81 (0)289-64-3380
FAX : +81 (0)289-62-5636

NSK Europe GmbH
www.nsk-europe.de
Elly-Beinhorn-Strasse 8
65760 Eschborn,
Germany
TEL : +49 (0)6196 77606-0
FAX : +49 (0)6196 77606-29

NSK United Kingdom Ltd
www.nsk-uk.com
Office 5, Gateway1000,
Arlington Business Park, Whittle Way,
Stevenage, SG1 2FP, UK
TEL : +44 (0)1438 310670
FAX : +44 (0)1438 310671


NSK Oceania Pty Ltd
www.nsk-inc.com
Unit 22, 198-222 Young St.
Waterloo, Sydney,
NSW 2017, Australia
TEL : +61 (0)2 8306 3000
FAX : +61 (0)2 8306 3030

NSK America Corp
www.nsk-inc.com
700 Cooper Court
Schaumburg, IL 60173,
USA
TEL : +1 847 843 7664
FAX : +1 847 843 7622

NSK Middle East
www.nsk-inc.com
Room 6EA-701, 7th Floor, East Wing No.6
Dubai Airport Free Zone,
PO Box 54316 Dubai, UAE
TEL : +971 (0)4-204-5507
FAX : +971 (0)4-204-5508

Specifications are subject to change without notice.

REF PR-D356E Ver.2 09.03.00.N



EX Series Handpieces
ISO E-Type External Water Spray

1:1 Direct Drive



MODEL **FPB-EC** ORDER CODE **Y110-158**

- For CA burs (ø2.35)
- Clean Head system
- Push Button Chuck



MODEL **NM-EC** ORDER CODE **Y110-161**

- Miniature head
- Latch type
- For CA burs (ø2.35 Short Shank)
- Max. 30,000 min⁻¹



MODEL **NAC-EC** ORDER CODE **Y110-148**

- For CA burs (ø2.35)
- Latch type
- Max. 25,000 min⁻¹



MODEL **FFB-EC** ORDER CODE **Y110-162**

- For FG burs (ø1.6)
- Clean Head system
- Ultra Push Chuck



MODEL **BB-EC** ORDER CODE **Y110-159**

- For CA burs (ø2.35)
- Latch type
- Max. 40,000 min⁻¹



MODEL **NRS2-EC** ORDER CODE **Y110-163**

- For CA burs (ø2.35)
- Clean Head system
- Push Button Chuck



MODEL **NBBW-EC** ORDER CODE **Y110-160**

- For CA burs (ø2.35)
- Latch type
- With Water Nozzle



MODEL **Spray Nozzle** ORDER CODE **C032-701**

- For NAC, BB, FPB, FFB (Optional)



MODEL **Spray Nozzle** ORDER CODE **C039-701**

- For NM (Optional)

135°C All handpieces and air motors can be sterilised up to 135°C.

EX Series Handpieces
ISO E-Type External Water Spray

1:5 Speed Increasing Contra-angles



MODEL **Ti95EX** ORDER CODE **C487**

- Solid Titanium Body
- For FG burs
- Max. 200,000 min⁻¹
- Ceramic Bearings
- Clean Head system
- Ultra Push Chuck

Prophylaxis Contra-angles

1:1 Direct Drive



MODEL **AR-EC (S)** ORDER CODE **Y110-164**

- For screw-in cups and brushes

4:1 Reduction



MODEL **AR-ER4 (S)** ORDER CODE **Y110-140**

- For screw-in cups and brushes



MODEL **AR-EC (K)** ORDER CODE **Y110-165**

- For snap-on cups



MODEL **EVA-ER4** ORDER CODE **Y110-141**

- For Eva tips
- Vertical Reciprocating Action
- Stroke 1.4 mm

135°C All handpieces and air motors can be sterilised up to 135°C.

Ti95EX can be washed in the thermodisinfectior.

EX Series Endodontic Handpieces
ISO E-Type External Water Spray

64:1 Reduction



MODEL **MP-ER64** ORDER CODE **Y110-142**
• Miniature head • For engine files
• 360° Rotation • Push Button Chuck



MODEL **MPA-ER64** ORDER CODE **Y110-144**
• Connection for Apex locator (Long shank files)
• Miniature head • For engine files
• 360° Rotation • Push Button Chuck



MODEL **MPAS-ER64** ORDER CODE **Y110-146**
• Connection for Apex locator (Short shank files)
• Miniature head • For engine files
• 360° Rotation • Push Button Chuck

100:1 Reduction



MODEL **NRS2-ER10** ORDER CODE **Y110-143**
• For engine files • 360° Rotation
• Clean Head system • Push Button Chuck
• 10:1 reduction handpiece & 10:1 reduction head



MODEL **TEP-ER10** ORDER CODE **Y110-145**
• For hand files • 90° Reciprocating Action
• Push Button Chuck



MODEL **Spray Nozzle** ORDER CODE **C032-701**
• For TEP (Optional)



135°C
All handpieces and air motors can be sterilised up to 135°C.

EX Series Straight Handpieces
ISO E-Type External Water Spray

1:1 Direct Drive



MODEL **EX-6B** ORDER CODE **H260-001**
MODEL **EX-6BM** ORDER CODE **H260-002**
• For HP burs (ø2.35) • Clean Head system



MODEL **EX-6D** ORDER CODE **H267-001**
MODEL **EX-6DM** ORDER CODE **H267-002**
• For HP burs (ø2.35) • Clean Head system


4:1 Reduction



MODEL **EX-5B** ORDER CODE **H262**
• For HP burs (ø2.35) • Clean Head system



Ergonomic Triangular Grip
EX straight handpieces have the ergonomic triangular cross section which provides a comfortable and relaxing hold for prolonged clinical procedures.



MODEL **Spray Nozzle** ORDER CODE **H019-711**
• For EX-6B, EX-5B (Optional)

EX-203C Air Motor Pack (EX-203C+EX-6B+NAC-EC)
ISO E-Type External Water Spray

1:1 Direct Drive



MODEL **EX-203C M4 Pack** ORDER CODE **Y160-020**
• Midwest 4 hole

MODEL **EX-203C B2 Pack** ORDER CODE **Y160-021**
• Borden 2 hole

- Contents
- Air Motor (EX-203C)
 - Straight handpiece (EX-6B)
 - Contra-angle handpiece (NAC-EC)

135°C
All handpieces and air motors can be sterilised up to 135°C.

EX-6BM and EX-6DM can be washed in the thermidisinfector.