



## Tri Auto mini

The new compact endodontic motor

### Advantages of Tri Auto mini at a glance

- New and innovative functions
- Rechargeable battery operation and cable-free
- Flexible speed range from 50 – 1,000 rpm
- 6 storage locations for treatment parameters
- Modularly expandable with Root ZX mini
- For all standard NiTi files

Thinking ahead. Focused on life.



**MORITA**

# Tri Auto mini – Convincing functions

The high-performance Tri Auto mini endodontic motor excels with its extremely high level of flexibility. With its compact dimensions and low weight, it fits perfectly into the hand and opens up a new degree of freedom in treatment possibilities due to its rechargeable battery operation. The user can adjust file speed, torque values and various safety functions at the touch of a button and in this way increase safety for patients. And the various parameters for individual indications can be stored for even easier handling.

In combination with the Root ZX mini apex locator, you have an endometric and treatment instrument of the very highest quality.



LC display colour selection



Special display setting for left-handers

## Safety

- Innovative functions
  - Auto Torque Reverse
  - Auto Torque Slow Down: Reduction of the rotational speed when approaching torque limitation
  - Visual indication of file position on display\*
  - Change of display colour when approaching the apex\*

\*) in combination with Root ZX mini

## Modular System

- Upgrade with Root ZX mini apex locator per data transfer cable
- Additional functions for optimum preparation
  - Length Measurement
  - Auto Apical Stop
  - Auto Apical Reverse
  - Auto Start / Stop
  - Auto Apical Slow Down

## Tri Auto mini - Technical Data

### Handpiece

Display: LCD  
 Weight: 78 g  
 Motor: 50 – 1,000 rpm  
 Torque: 3.9 Ncm  
 Battery (rechargeable): 3.7 V lithium-ion  
 Diameter: 28 mm / Height: 150 mm

## Flexibility

- Tri Auto mini is suitable for all standard NiTi files
- Rotational handpiece (250°) for optimum visibility in all quadrants
- Speed range from 50 – 1,000 rpm
- Suitable for both right-handers and left-handers
- Battery-operated (rechargeable) and cable-free

## Treatment Comfort

- Easy-to-handle angle piece with extremely small head (9mm) for high patient comfort
- Clear, energy-saving LC coloured display
- Selectable signal tones
- Variable file speeds
- Variable torque values
- 6 storage locations for treatment parameters



Distributed by  
 J. Morita Europe GmbH  
 Justus-von-Liebig-Str. 27a, 63128 Dietzenbach, Deutschland  
 Tel.: +49 (0)60 74 / 836 - 0, Fax: +49 (0)60 74 / 836 - 299  
 www.JMoritaEurope.com

Manufactured by  
 J. Morita Mfg. Corp.  
 680 Higashihama Minami-cho, Fushimi-ku, Kyoto, 612 - 8533 Japan  
 Tel.: +81 (0)75 / 611 - 2141, Fax: +81 (0)75 / 622 - 4595  
 www.JMorita.com

